

PARAFFIN WAX BATH

Size: 26" x 18" x 16", Tank Size (18"x13"x8")
 Size: 23" x 13" x 14", Tank Size (14"x9"x7")
 PARAFFIN WAX

LCS 240
 LCS 241
 LCS 251

Internal wax tank made of SS sheet. The outer tank & top is made of MS sheet and the edges of the bath has a laminated Wooden rim, provided with special wax immersion type Heater & thermostat to control temperature, comes up with 5 kg. Paraffin wax.



HYDROCOLLATOR 4 PACKS Heavy Duty with 4 Packs.

LCS 262

HYDROCOLLATOR 8 PACKS Heavy Duty with 8 Packs.

LCS 271



LCS 262
 Made of stainless steel (Size 48cm x 38cm x 59cm), (Tank Size 35cm x 27cm x 29cm)
 Double walled with 1000 watts heater & thermostat for heat control comes up with 4 steam packs.

LCS 271
 Made of stainless steel (Size 63cm x 43cm x 79cm), (Tank Size 46cm x 31cm x 33cm),
 Double walled fitted with 2000 watts heater, and thermostat for heat control comes up with 8 steam packs

HYDROCOLLATOR PACKS: Available in 4 sizes.

- A-Large : size 11"x21"
- B-Medium : size 10"x12"
- C-Small : Size 07"x09"
- D-Cervical : size 09"x26"

LCS 281

Pack contains silica crystals having good amount of heat retention.

