

NEW STYLE

AYU-250W-ANL

**ELECTRO SURGICAL UNIT
(250W)**

- Light Weight
- Elegant Design
- Simple Operation
- The Reliable Economical



Ayu-250W-ANL

Electro Surgical

NEW STYLE

AYU-250W-ANL

ELECTRO SURGICAL UNIT

(250W)

www.lcsindia.com

HIGH FREQUENCY ELECTRO SURGICAL UNIT

microprocessor control, innovative design, together with high reliability makes 250W High Frequency Surgical Units the most favorable choice for doctors.

In the used surgical diathermy, the technical ability of the unit combined with plain safety are of prime importance-both for the patient and the user. There must be no possibility of accidental, undesirable results in the use of Ayunix Electro Surgical Units been developed.

The Ayu 250W-ANL are newly developed Electro surgical Units of particular safety and reliability. The Ayu 250W as the most refined unit of its kind in the world.

SAFETY IN USE


The technical conceptions of the Ayu 250W-ANL meet the requirements of highly refined, Electro Surgical techniques all national and international safety standards are conformed to or exceeded. It is not possible for the unit to deliver a higher output than has been set on the control. The micro function allows for an instant reduction in the output for more delicate application such as neurological or micro surgery.

OUTPUT CHARACTERISTICS

	MODE	MAXIMUM P-P VOLTAGE	RATED LOAD (OHMS)	MAXIMUM POWER (WATTS)	CREST FACTOR* (TYPICAL)
Monopolar Cut	Pure	1400	500	250	1.5
	Blend	1800	300	180	2.5
Monopolar Coag.	Force	2200	300	180	5.0
	Soft	1800	300	150	4.5
Bipolar Cut	Pure	500	100	100	1.5
	Blend	700	100	70	2.0
Bipolar Coag	Force	735	100	70	2.5
	Soft	600	100	50	2.0

Physical Characteristics - L x W x H - 430x405x280 mm
Net Weight - 5Kg.200Grm. Gross Weight - 10Kg. 500Grm.

Standard Accessory



Patient Plate

Hand Switch

Unipolar Handle

Knife Electrode (Angled)

Knife Electrode (Straight)

Needle Electrode

Ball Electrode 6mm

Loop Electrode 5mm

Bipolar Forceps

Unipolar Footswitch



Ayu-250W-ANL
Electro Surgical

LIFE CARE SYSTEMS

K-708, Sector-23, Sanjay Nagar,
Ghaziabad-201002 Delhi NCR
Cell : +91 9312247417
Cell : +91 9654875195
E-mail : sales@lcsindia.com
Website : www.lcsindia.com