

ADVANCE

AYU-43AB

ELECTRO SURGICAL UNIT (400W)

- Light Weight
- Elegant Design
- Simple Operation
- The Reliable Economical

www.lcsindia.com



High Frequency Electro Surgical Unit With Separate Boiler Footswitch

microprocessor control, innovative design, together with high reliability makes 400W High frequency Surgical Units the most favorable choice for doctors.

In the used surgical diathermy, the technical ability of the unit combined with plain safety are of prime importance-both for the patient and the user; There must be no possibility of accidental, undesirable results in the use of Ayunix Electro Surgical Units been developed.

The Ayu 400W-ANL are newly developed Electro surgical Units of particular safety and reliability. The Ayu 400W as the most refined unit of its kind in the world.

SAFETY IN USE


The Technical conceptions of the Ayu 400W-ANL meet the requirements of highly refined, Electro Surgical techniques all national and international safety standards are conformed to or exceeded. It is not possible for the unit to deliver a higher output than has been set on the control. The micro function allows for an instant reduction in the output for more delicate application such as neurological or micro surgery.

OUTPUT CHARACTERISTICS

	MODE	RETED LOAD (OHMS)	MAXIMUM POWER (WATTS)
Monopolar Cut	PURE	300	400
	BLEND-I	300	300
	BLEND-II	300	200
Monopolar Coag.	DESSICCATE	300	150
	FULGURATE	300	250
	SPRAY	300	120
Bipolar Cut	MACRO	100	100
	MICRO	100	70
Bipolar Coag	MACRO	100	70
	MICRO	100	50

Physical Characteristics - L x W x H - 430x405x280 mm
Net Weight - 5Kg 500Grm. Gross Weight - 10Kg. 900Grm.

Standard Accessory



LIFE CARE SYSTEMS

K-708, Sector-23, Sanjay Nagar,
Ghaziabad-201002 Delhi NCR
Cell : +91 9312247417
Cell : +91 9654875195
E-mail : sales@lcsindia.com
Website : www.lcsindia.com