

CONTINUOUS PASSIVE MOTION**LCS 231**

LCS CPM is designed for lower limb. It is based on solid-state circuitry, silent, smooth & powerful servo motor for continuous operation, instant start/stop.

APPLICATION

Continuous Passive Motion therapy immediately used after surgery can give several important benefits to patient like maintenance or improvement in the range of motion, reduction in pain and discomfort.

SPECIFICATION: -

Dimension (L*H*D)	1100mm, 224mm, 185mm
Weight	16 Kg. (Approx)
Operating Voltage	220v 50Hz
Absorption	60-Watts Max.
Fuses	1 Amp
Room Temperature	10° to 40°c
Moisture	10% to 80%
Notation Range	10° to 120°
Digital Timer	01 ~ 99 Minutes Adjustable
Patient Safety	Available

**TRACTION TABLE****LCS 302**

Fixed Height 4 section Traction Table for horizontal Cervical & Lumbar Traction Treatment, comes up with adjustable stool. We offer smooth and comfortable traction for patients.

SPECIFICATION: -

Height	27 Inches (68.58 cm)
Width	24 Inches (60.96 cm)
Length	94 Inches (238.76 cm)

**FEATURES INCLUDE**

- Friction-free lumbar section.
- Comes up with adjustable stool.
- Cushion with comfortable slow recovery foam.
- Lever system to allow front OR back seat angle adjustment.
- One year limited warranty.

